

崩壞試驗器 Disintegration Tester



特點：

- 本機可試驗藥品於定溫下之崩散程度，其振盪速度採無段變速方式。反應水槽內的溫度可依需要任意設定並採自動恆溫循環裝置，並加裝有漏電斷電防護裝置及過溫防護裝置，以確保使用安全。懸籃採SUS 316製造。

Feature :

- It is used to test the disintegration condition of chemicals under the set temperature. It's shaking speed is non-limitary variation. The temperature of water baths can be set to a need one, in special, it can keep up the set temperature automatically and the circuit breaker device can prevent from leakage current. Thus, it can provide safety when use. The basket are made of SUS316.

Model	DG-2B
懸籃 Basket	2
懸籃動作 Movement of Basket	Amplitude 55mm, Vertical agitation 5 ~ 40 times/min
驅動馬達 Driving Motor	6W
加熱器 Heater	600W
尺寸 Dimension (mm)	470 x 325 x 555
重量 Weight (kg)	22
電源 Power	AC 220V, 50/60Hz

配件 / Accessories :



懸籃
Basket
2 PCS



燒杯
Thick Beaker
2 PCS



輔助筒
Granual Cylinder
12 PCS



輔助盤
Disk
12 PCS