

敝公司提供電熱包系列產品，被廣泛使用於實驗室及產業界，而其表面最高溫度達攝氏450度，並接受客戶特殊規格品製作。

### Mantle Heater For Flask 標準型電熱包



MODELS	CAPACITY	POWER
MN-50	50 ml	60 w
MN-100	100 ml	80 w
MN-250	250 ml	200 w
MN-500	500 ml	300 w
MN-1000	1000 ml	400 w
MN-2000	2000 ml	500 w
MN-3000	3000 ml	600 w
MN-5000	5000 ml	700 w
MN-10L	10 L	1300 w
MN-12L	12 L	1300 w

### Mantle Heater For Flask With Controller 附調溫器電熱包



MODELS	CAPACITY	POWER
MN-50C	50 ml	60 w
MN-100C	100 ml	80 w
MN-250C	250 ml	200 w
MN-500C	500 ml	300 w
MN-1000C	1000 ml	400 w
MN-2000C	2000 ml	500 w
MN-3000C	3000 ml	600 w
MN-5000C	5000 ml	700 w
MN-10LC	10 L	1300 w
MN-12LC	12 L	1300 w
MN-20LC	20 L	2400 w
MN-22LC	22 L	2400 w
MN-50LC	50 L	3200 w

### Heating Tapes 電熱帶



MODELS	WIDTH	LENGTH	POWER
TIP-2	25 mm	2 ft	100 w
TIP-4	25 mm	4 ft	200 w
TIP-6	25 mm	6 ft	300 w
TIP-8	25 mm	8 ft	400 w
TIP-12	25 mm	12 ft	600 w
TIP-16	25 mm	16 ft	800 w

電熱帶長度單位為英尺；最高表面溫度400°C

+°C

Heating Mantles