

# PORTABLE VISCOMETER B-ONE



This viscometer  
follow you  
everywhere



## FEATURES AND BENEFITS

- Autonomy more than one hour.
- Charge indicator.
- Direct measure with time to stop.
- Usable with standard or rack stand (option).
- User and locked mode.
- Data recording and USB transfer.
- Torque gage on display.
- Display of viscosity limits according to mobile and speed.
- Speed or shear rate control.

## AVAILABLE INSTRUMENT

Reference	Designation	Viscosity range (mPa.s)
T600700	PORTABLE B-ONE VISCOMETER WITHOUT MEASURING SYSTEM	200(15**) to 240M*

M for million / \* With MS-RV (RV2 to RV7). Not included. / \*\* Can reach 15 mPa.s with MS BV.

WHAT'S INCLUDED WITH DEVICE ? (according to part number)

1 Carrying case / 1 User manual / 1 Calibration and checking certificate / 1 Harness.

## COMPATIBLES MEASURING SYSTEMS



MS-RV, MS-LV, MS-BV, MS-VANES,  
MS-KREBS, MS-SV, MS-ULV, MS-DIN.

## COMPATIBLES TEMPERATURE CONTROLS



EVA LR, EVA MS DIN,  
RT-1 PLUS.

## OPTIONS & ACCESSORIES



- Standard or rack stand.
- HELIPRO stand.
- Small volume package.
- THERMOCELL package.
- Bayonet spindle rack.
- Barcode reader.
- Touch screen protector.
- Viscosity standards oil.

## SPECIFICATIONS



### Type of instrument

Rotating springless viscometer with 7" Touch screen

### Rotation speeds

Unlimited number of speeds between 0.3 and 250 rpm

### Torque range

From 0.05 to 13 mNm

### Accuracy

+/- 1 % of the full scale

### Repeatability

+/- 0,2 %

### Display

Viscosity (cP / Poises or mPa.s / Pa.s)  
Speed – Shear rate – Torque (% - mNm) –  
Shear stress – Time

### Languages

French/English/Russian/Spanish/Turkish/  
German/Italian/Korean/Portuguese/Ukrainian

### Supply voltage

90-240 VAC 50/60 Hz

### Dimensions et poids

Head : Ø 90 mm  
Box D60 x H260 x W120  
Weight: 2 kg

