

## 压片机 PP 25

### 基本信息

稳定而有效的片剂是进行可靠X射线荧光分析的重要前提。本款液压手动式压片机PP25是一台设计紧凑实验室设备，操作起来相当简单而且安全。它能产生的最大压力为25吨，尤其适合于制备用于X射线荧光分析的固体样品。其制压而成的样品片同时具备了高质量和高稳定性。压力塞的压力可通过测压计清楚地读出。PP25的可拆式压模有多种直径型号提供。



### 应用案例

原料, 水泥, 炉渣, 矿物, 矿石, ...

### 产品优势

- 生产高质量的、稳定的颗粒
- 操作简单、安全
- 提供多种尺寸的压片模具
- Can be evacuated
- 桌面型仪器，设计紧凑

### 特点

应用	为光谱分析压片
应用领域	化学 / 合成材料, 地质 / 冶金, 建筑原料, 环境 / 资源回收利用, 玻璃 / 陶瓷
最大压力	25 吨, 液压
压力	1 - 25 t (10 - 250 kN)
压片模具	∅ 32 / ∅ 40 mm
机体尺寸 (宽x高x纵深)	380 x 620 x 300 mm
净重	50公斤
标准	CE
请注意:	
受样品材料性质和仪器配置/设定的影响	

### 视频链接

<http://www.retsch.cn/pp25>

### 功能原理

样品被填充于两片经抛光处理的压力片之间。可活动的传压轴的一端突出于模具的外壁，由液压系统产生的压力经由此处传到样品。

### 订单数据

压片机PP25

## 压片机 PP 25

( 请单独订购压片工具 )

20.750.0006

压片机PP25，液压，手动型

### PP25克抽真空型压片工具

22.458.0016

Pressing tool for aluminum cups Ø 32 mm (also suitable for free pressing)

22.458.0017

Pressing tool for aluminum cups Ø 40 mm (also suitable for free pressing)

further pressing tools upon request

### PP25配件

22.005.0001

铝杯，倾斜壁，适用于直径32毫米的压片，1000件

22.005.0002

铝杯，倾斜壁，适用于直径40毫米的压片，600件

22.458.0006

Aluminum cups, straight walls, for pellets with 40 mm diameter, 1000 pieces

22.440.0001

Licowax® C micropowder, 250 g (not for steel rings)

22.440.0003

纤维素黏合剂Spektromelt® C20，1公斤