

Mini UV Viewing Cabinets

The mini Chromato-Vue® UV Viewing Cabinets are designed with a metal base and plastic top. The **C-15G** cabinet includes a UVGL-58 UV lamp. Lamps for other mini cabinets are ordered separately.

- The bottom panel of the cabinet is removable
- The lamp fits snugly on top of the cabinet
- A contrast control filter absorbs UV energy and protects the eyes
- Soft rubber viewport provides a comfortable viewing area with space for eyeglasses
- Lightweight access curtain blocks out external light from the cabinet interior

C-10P Mini UV Cabinet accommodates the 6-watt Rechargeable Lamps which are excellent for field use applications (see page 6 for Rechargeable Lamps).



Mini Cabinet Dimensions:
 Interior: 8.75L x 10W x 7H in. (222 x 254 x 178mm)
 Exterior: 9L x 10.5W x 12H in. (229 x 267 x 305mm)



C-15G Mini UV Cabinet Combines the C-10 Cabinet with UVGL-58 multi-band lamp with shortwave and longwave.



C-10E4 Mini UV Cabinet Designed for use with the 4-watt EL Series Lamps (see pages 4 - 5).

Ordering Information

C-10 Mini Cabinet for 4 and 6-watt PL Series UV Lamps	95-0072-01
C-10P Mini Cabinet for 6-watt Portable Rechargeable Lamps	95-0072-02
C-10E4 Mini Cabinet for 4-watt EL Lamps	95-0072-08
C-10E6 Mini Cabinet for 6-watt EL Lamps	95-0072-09
C-15G Mini Cabinet with UVGL-58 UV Lamp	95-0072-06 (115V) 95-0072-07 (230V)

Multifunctional UV Cabinet

C-65 Chromato-Vue Cabinet accommodates one or two EL Series 8-watt lamps ordered separately. Several lamps are listed below. For additional lamps, see pages 4 - 5. Use of two lamps illuminates the interior uniformly from two sides of the cabinet.

- The bottom panel of the cabinet is removable
- Includes an interior overhead white light
- Access curtain blocks ambient light

Ordering Information

C-65 Mid-Size UV Viewing Cabinet	95-0257-01 (115V) 95-0257-02 (230V)
---	--

C-65 Dimensions:
 Interior: 17L x 14W x 6H in. (533 x 330 x 400mm)
 Exterior: 19.75L x 16.5W x 15.75H in. (602 x 419 x 400mm)

EL Series 8-watt lamps can be used with the C-65 cabinet or removed for handheld use (see p. 4-5 for more EL Lamps).

- Lamps are made with a durable, scratch-resistant metal
- Lamps include a universal plug, which connects directly to the C-65 Cabinet, plus a supplementary power cord for use of the lamp separately from the cabinet



C-65 UV Cabinet is designed to accommodate one or two of the EL Series 8-watt Lamps

Ordering Information - EL Lamps

UVL-28 365nm/365nm UV	95-0248-01 (115V) 95-0248-02 (230V)
UVS-28 254nm/254nm UV	95-0249-01 (115V) 95-0249-02 (230V)
UVLS-28 (2UV) 365nm/254nm UV	95-0201-01 (115V) 95-0201-02 (230V)
UVLMS-38 (3UV) 365nm/302nm/254nm UV	95-0252-01 (115V) 95-0252-02 (230V)