



# VC 505

500 Watts

Processing Volume: 200µl – 1L\*



The **VC 505** is a cost effective alternative to the VCX model. This 500 watt system offers the same processing power and choice of accessories as the VCX500 without temperature monitoring or energy set point programming.

A standard ½" diameter probe with replaceable tip is included and many optional accessories are available to meet the needs of any application.

## FEATURES:

- Programmable
- Variable power output
- 10 hour timer
- Pulse mode
- Wattage and Energy display



(stand is optional)

Ordering information: Model no. **VC 505**  
500 watt ultrasonic processor with ½" diameter probe and tool kit.

## SPECIFICATIONS

Power Output	500 watt maximum	Converter Cable Length	Part #201-0300 6 ft. (1.8m)
Frequency	20kHz		
Dimensions	H x W x D: 9.25 x 7.5 x 13.5" (235 x 190 x 340mm)		
Ultrasonic Converter	Part #CV334	Tool Kit	Part #381-0005 Includes: 2 spanner wrenches and 1 open-end wrench
Diameter	2.5" (64mm)		
Length	7.25" (183mm)		
Standard Probe	Part #630-0220	Optional:	
Tip Diameter	0.5" (13mm) with replaceable tip	Stand with Clamp	Part #830-00459
Length	5.5" (139mm)	Sound Enclosure	Part #830-00427
Material	Ti6Al4V titanium alloy		

\*The actual processing volume of a specific probe or accessory is application specific.